

SEED/GRAIN DIVIDER

BOERNER TYPE THE PLT - 148(b) SERIES

APPLICATIONS

- Seed Grain Divider Boerner Type is mainly used for local grain inspection.
- Used for dividing seed sample into two equal halves, as per FCI and B.I.S Requirements.

SALIENT FEATURES & CONSTRUCTION

- This divider is made to meet **F.C.I. and B.I.S Requirements.**
- Especially for free flowing grain lots.
- Smooth finish without joint
- Brass Body which makes it rust resistant.
- Two Brass sample pans for easy transport of divided samples.
- Model with 36 openings or pockets
- Gravity operated.
- No Electric requirement & easy operation.



An ISO 9001: 2008 Certified Co.

