

WATER BATH RECTANGULAR

SINGLE WALLED **THE PLT - 107 SERIES**

APPLICATIONS

- Single Walled Water Bath generally used for evaporating liquids for boiling the material.
- Single Walled Water Bath is low cost options for heating / boiling ideal tool for Laboratories, Pharmaceutical, Food & Beverage Industries and Institutions.

SALIENT FEATURES & CONSTRUCTION

- Single walled made of heavy gauge stainless steel.
- Provided with concentric rings (Removable) of 75mm hole.
- Concentric rings can be removed to accommodate containers of various sizes.
- Heating by immersion heaters.
- Temperature range ambient up to 90°C.
- Operating Voltage: 220 Volts AC (50 Hz)

Model No.	WBR-1	WBR-2
Internal dimensions (L X W X D)	300 x 250 x100mm	350 x 400 x100mm
Capacity	8 Ltr	15 Ltr
No. of Holes	6	12
Holes Dia.	75 mm	75 mm
Heating Load	1.0 KW	2.0 KW

