

WATER STILL AUTOMATIC

IN ALL GLASS

THE PLT - 158 SERIES

APPLICATIONS

- These water distillation plants fully automatic have been designed to meet with the requirements of best quality of distilled water. These are best suitable for removing all types of impurities Organic, Inorganic, Gaseous, Bacterial and Pyrogenic leaving a final pure distillate. (As per NPL Patent.)

SALIENT FEATURES & CONSTRUCTION

- The entire range is fabricated using premium quality glass material (**BOROSIL GLASS**).
- Produces Ultra pure, Pyrogen free distilled water.
- Fitted with Automatic low level water cut-off device.
- The distillation apparatus consists of flask with heating elements embedded in glass and fused in spiral type coil internally of the bottom and tapered round glass, joints at the top double walled condenser with B-40/B-50 ground glass joints.
- Constant level arrangement for water. Individual pilot lights on the front panel to show the status of water level and water flow.
- Complete with metallic stand rings and clamps to hold flasks condenser electric board with switch, Neon light indicator and connecting cords.
- Supplied with cord and plug to work on 220/230 volts A.C. supply.



Model No.	WSG-3	WSG-5	WSG-10
Flask Capacity	3 ltr	5 ltr	10 ltr
Available Models	Single Stage, Double Stage, Triple Stage		
Operating Voltage	220 Volts AC		

Tanco[®]

An ISO 9001: 2008 Certified Co.

